English

1. APPLY THE SOCK PROVIDED

2. PREPARE THE BOOT

Remove the foam from the boot (Charcot Conformer: Also remove the black plastic PTB plate.)

3. APPLY THE BOOTIE

(A) Slip the foot into the bootie with the toes beneath the tongue. Slide the heel of the foot all the way to the back of the bootie. Make certain the foot is centered from the left to right on the white innersole section.

(B) Close the foam on the forefoot.

4. SECURE THE BOOTIE

(A) Wrap the foam section around the leg by tucking in the long end of the foam cuff.

(B) Adjust the tension on the wrap, and close the hook and loop closure.

(C) If the cuff is too long, trim with scissors.

5. APPLY THE BOOT

Make sure the plastic covers are still on the boot uprights. Slip the bootie into the boot, position it evenly front to back. Stand and place full body weight on the boot for 30 seconds. NOTE: When properly positioned, the white innersole should fit completely inside the black midsole.

6. SECURE THE UPRIGHTS

(A) While still bearing weight and with the foot pointed directly forward, align the uprights with the sides of the leg parallel to the leg bones.

When the uprights are properly aligned with the leg, remove the plastic from the uprights and press the uprights against the side of the leg to engage the hook material to the bootie.


(B) Save the plastic covers for reapplication after washing the bootie.

7. SECURE THE FOREFOOT STRAPS

While standing, wrap the forefoot straps over the top of the foot, through the D-rings and pull back to tension. If necessary, trim the excess strap length with scissors.

NOTICE TO USERS

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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(C) Repet steps A and B for the remaining straps.

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